

## OPERASIONAL DAN MANUAL CB-001-D

No	TYPE	KARAKTERISTIK	SPESIFIKASI KINERJA	ANALISA RESIKO ALAT
1	CB-001 D	Sistem pengaturan ketinggian bed menggunakan lock pin dan standar level ketinggian diatur dalam 3 tingkat ketinggian (knock down),	Pengoperasian secara manual dengan menaikkan dan menurunkan bed sesuai posisi lock pin	Tidak ada
		Cara penyesuaian ketinggian dilakukan dengan manual (knock down) dengan bantuan tenaga manusia.	Penyesuaian ketinggian dilakukan dengan secara manual dengan tenaga manusia saat kondisi bed kosong (tanpa beban)	
		Sistem pengaturan level ketinggian dengan menggunakan lock pin yang dirancang khusus yang berfungsi sebagai stopper (kaitan), untuk menahan beban saat bed digunakan	Lock dirancang untuk menahan beban pasien sesuai standar yang ditetapkan	
		Menggunakan penggerak roda caster dengan double lock system untuk keamanan saat parkir	Double lock terdiri dari lock caster kanan depan dan kiri belakang	Pemakaian yang terus menerus dan kasar akan membuat sistem lock cepat rusak

### CB-001 D (WITH CASTERS)



#### CASTER OPERATION SYSTEM

Independent lock caster (one at front and one at back of the bed) provides stability and secure parking.



#### PRODUCT SPECIFICATIONS

**Operating System** Three adjustable Hi-Lo (knock down) / **Accessories** Side Rail (CSR-001, CSR-002, CSR 003)

**Main Frame Material** Steel / **Head & Foot Frame Material** Steel / **Weight (kg)** Approx. 60 / **Frame Finishing** Powder Coating

**Frame Colour** Ivory / **Dimension W x D x H (mm)** 930 x 2129 x 932

#### PRODUCT DIMENSION

