

CB 3003 PLUS



Standard Function

- Operating System : 3 Crank operation for High-Low, Head Rest and Foot Rest positions
- Position : - High - low: Manual operation with Crank
- Head rest: Manual operation with Crank, tilting from 0 to 80 degrees
- Foot Rest: Manual operation with Crank, tilting from 0 to 80 degrees
- Primary : - Simple design, no sharp corners (smooth surface)
- Easy operating mechanism, with Crank System
- Sturdy and anti-rust
- Wheels that are easy to move and position, with lock - unlock system
- Lock System : Double Lock System on the Wheels
- Accessories : IV Pole (2 Type Option) and Side Rail / Folding Side rail
Bed Side Cabinet, Over Bed Table, and Foot Step
Aimate Cpro Mattress / Urethane



**PERHATIAN !
PENGOPERASIAN HANDLE CRANK**

- Sebelum memutar handle crank, harap perhatikan :
Tangan Anak-anak, Hewan peliharaan, alat-alat dan bahan-bahan
Befores using the crank handle, please watch :
Hands, children, pets, other people, materials.
- Setelah memutar Handle, pasang kembali pada posisi terkunci untuk
menghindari jatuh Dae-tiba
After turning the crank handle, put back on locked position in order to avoid sudden falling down.
- Untuk menghindari terkane kaki, pasang Handle Crank pada posisi awal
(terkunci). Setelah selesai, maka Bed akan jatuh
In order to moving the, put the crank handle back on first position (locked). It is to, the bed weight fall
down.

Caution Handle operation



IZIN EDAR KEMENKES

PERHATIAN !

- Jangan menempatkan tangan
dibawah tempat tidur.
Do not place your hand under the bed
- Tempat tidur hanya digunakan
untuk satu orang.
Capacity for one person only

CAUTION LABEL

Lot No :
Type :

PT. Chitose Internasional Tbk
Jln. Industri III No.5 Utama Cimahi
40533, Jawa Barat Indonesia

Showroom - Marketing Office
Distributor Center Jln. H.M.S.
Mintaredja Baros Cimahi 40521,
Indonesia Ph. +62 22 206 4777

Manufactured by Chitose

MANUFACTURE LABEL AND LOT NUMBER

- 1 Lot Number and Product Type
- 2 Manufacture Information
- 3 Safe Workload Sign and Maksimum Weight Sign
- 4 Use Appropriate mattress Sign
See Manual Instruction Sign
- 5 Applied Part Sign
Do not Dispose as Household Waste Sign
Indoor Using Only Sign