

CHITOSE Hospital Bed [trademark]

CB3003 plus



FEATURES & PROCEDURES



Head/foot board panel dismantle

Head/foot board panel installation

Irrigator bar installation (applicable for every nursing bed)



An unique eye catching foldable crank is designed to give high-low bed adjustable

Independent lock caster (one at front and one at back of the bed) provides stability and secure parking

Fungsi Standar

Sistem Operasi : 3 Operasi Crank untuk posisi Tinggi-Rendah, Sandaran Kepala & Sandaran kaki

Posisi : - Tinggi - Renda : Digerakkan manual dengan Crank
- Sandaran Kepala : Digerakkan manual dengan Crank, mulai dari 0° sampai 80°
- Sandaran Kaki : Digerakkan manual dengan cank, mulai dari 0° sampai 80°

Keutamaan : - design sederhana, sudut tidak tajam (permukaan halus)
- Mekanisme pengoperasian sederhana, menggunakan sistem crank
- Kokoh dan anti karat
- Roda yang mudah digerakkan dan ditempatkan, dengan sistem kunci - tanpa kunci

Sistem Penguncian : Sistem penguncian ganda pada roda caster

Assesories : tiang IV (2 tipe pilihan) dan Side Rail biasa/ Side rail lipat, lemari samping, Meja atas kasur, Pijakan kaki, matras Cpro/ Busa Urethane

Ukuran Produk

Ukuran

P x L x Tmin-max : 2230 x 1030 x 760~1005

Berat Produk

Berat : ± 90KG

Kapasitas beban : ± 400KG (Statis)
& ±170KG (Dinamis)

Ukuran Kemasan : 2050 x 1070 mm x 190

Lot No :
Type :

PT. Chitose Internasional Tbk
Jln. Industri III No.5 Utama Cimahi
40533, Jawa Barat Indonesia

Showroom - Marketing Office
Distributor Center Jln. H.M.S.
Mintaredja Baros Cimahi 40521,
Indonesia Ph. +62 22 206 4777

Manufactured by Chitose

MANUFACTURE LABEL AND LOT NUMBER

Chitose®
MANUFACTURER LABEL

SUCOFINDO SNI
PCS 01095.01 SNI IEC 60601-2-52:2014
LSPR-022-IDN

SNI LABEL

Chitose bed
KEMENKES RI AKD:

IZIN EDAR KEMENKES

bahan dasar : sticker warna putih
posisi penempatan pada Frame

