


CB 3012 PLUS



NURSING BED
Healthy Therapy Product



Standard Function
 Operating System : 2 Crank operation for High-Low and Head Rest positions
 Position : - High – low: Manual operation with Crank
 - Head rest: Manual operation with Crank with tilting from 0 to 80 degrees
 Primary : - Simple design, no sharp corners (smooth surface)
 - Easy operating mechanism, with Crank System
 - Sturdy and anti-rust
 - Wheels that are easy to move and position, with lock - unlock system
 Lock System : Double Lock System on the Wheels
 Accessories : IV Pole (2 Type Option) and Side Rail / Folding Side rail
 Bed Side Cabinet, Over Bed Table, and Foot Step
 Aimate Cpro Mattress / Urethane



CRANK LABEL



MANUFACTURER LABEL




Caution Handle operation



IZIN EDAR KEMENKES



CAUTION LABEL



MANUFACTURE LABEL AND LOT NUMBER

- 1 Lot Number and Product Type
- 2 Manufacture Information
- 3 Safe Workload Sign and Maksimum Weight Sign
- 4 Use Appropriate mattress Sign
See Manual Instruction Sign
- 5 Applied Part Sign
Do not Dispose as Household Waste Sign
Indoor Using Only Sign