

CB 3003 D

OPTIONS



CB-3003 D ST

FEATURES & PROCEDURES



Head/foot board panel dismantle

Head/foot board panel installation

OPERATION SYSTEM

A unique eye catching foldable crank is designed to give high-low bed adjustment.



Crank operation system

Side rail compartment (available in every nursing bed)

CASTER OPERATION SYSTEM

Independent lock caster (one at front and one at back of the bed) provides stability and secure parking.



PRODUCT SPECIFICATIONS

Operating System Adjustable Hi-Lo, back raise, and foot raise manual by crank, adjustable back raise with cylinder gas control

Accessories Side Rail (CSR-001, CSR-002, CSR 003) & Folding Side Rail (FSR-001)


Lock System Double Lock System (CB-3012 D, CB-7012 D), Four Lock System (CB-3012 D4)

Weight (kg) Approx. 90 / **Frame Finishing** Powder Coating

Frame Colour Ivory / **Frame Material** Steel

Dimension (mm) Bed (W 932, L 2141, H 654-994, H/F Board (W 932)

PRODUCT DIMENSION



Logo Perusahaan




Keterangan :
1. Bahan dasar sticker warna putih
2. Posisi penempelan pada frame

Label AKD




Label SNI



Alamat pembuat produk

Lot No & Tipe

Beban maksimum



Beban maksimum :

Beban Statis 245 kg

Beban Dinamis 175 kg

PT Chitose Internasional, Tbk		PT Indomedik Niaga Perkasa
Factory and Head Office	Showroom- Marketing Office	Distributor Head Office
Jln. Industri III No.5 Utama Cimahi 40533, Indonesia Ph.+62 22 603 1900 Fax.+62 22 603 1855	Jln. H.M.S. Mintaredja Baros Cimahi 40521, Indonesia Ph. +62 222 206 4777	Kem Tower Jl. Landas Pacu Barat No.10, RT 02 RW 10 Gn Sahari Sel Kemayoran Jakarta, DKI 10610 Ph.+62 21 657 04122 Fax.+62 21 657 04121

