

CHITOSE Hospital Bed [trademark] CB3003 plus



FEATURES & PROCEDURES



Head/foot board panel dismantle

Head/foot board panel installation

Irrigator bar installation (applicable for every nursing bed)



An unique eye catching foldable crank is designed to give high-low bed adjustable

Independent lock caster (one at front and one at back of the bed) provides stability and secure parking

Standard Function

Operating System : 3 Crank operation for High-Low, Head Rest & Foot Rest positions

Position : - High – low: Manual operation with Crank
 - Head rest: Manual operation with Crank, tilting from 0 to 80*
 - Foot rest : manual operation with Crank, tilting form 0* to 80*

Primary : - Simple design, no sharp corners (smooth surface)
 - Easy operating mechanism, with Crank System
 - Sturdy and anti-rust
 - Wheels that are easy to move and position, with lock - unlock system

Lock System : Double Lock System on the Wheels

Accessories : IV Pole (2 Type Option) and Side Rail / Folding Side rail
 Bed Side Cabinet, Over Bed Table, and Foot Step
 Airmate Cpro Mattress / Urethane

Product Dimension

Dimension

L x W x D : 2230 x 1030 x 760~1005


Product Weight

Weight : ± 90KG

Load Capacity : ± 400KG (Static)

& ±170KG (Dynamic)

Packing Dimension : 2050 x 1070 mm x 190




Lot No :

Type :

PT. Chitose Internasional Tbk
 Jln. Industri III No.5 Utama Cimahi
 40533, Jawa Barat Indonesia

Showroom - Marketing Office
 Distributor Center Jln. H.M.S.
 Mintarejda Baros Cimahi 40521,
 Indonesia Ph. +62 22 206 4777

Manufactured by Chitose



Chitose[®]

MANUFACTURER
 LABEL



PCS 01095.01 SNI IEC 60601-2-52-2014
 LSP-022-IDN

SNI LABEL

MANUFACTURE LABEL AND LOT NUMBER

Chitose bed

KEMENKES RI AKD: _____

IZIN EDAR KEMENKES

bahan dasar : stiker warna putih
 posisi penempatan pada Frame