

CB3012 plus



FEATURES & PROCEDURES



Head-foot board panel dismantle

Head-foot board panel installation

Trigero can be tilted and adjustable for every nursing care



An unique eye-catching foldable crank is designed to give high-low bed adjustable.



Embossment lock (upper part of foot and end of back of the bed) provides stability and secure parking.

Standard Function

- Operating System : 2 Crank operation for High-Low and Head Rest positions
- Position : - High – low: Manual operation with Crank
- Head rest: Manual operation with Crank, tilting from 0° to 80°
- Primary : - Simple design, no sharp corners (smooth surface)
- Easy operating mechanism, with Crank System
- Sturdy and anti-rust
- Wheels that are easy to move and position, with lock - unlock system
- Lock System : Double Lock System on the Wheels
- Accessories : IV Pole (2 Type Option) and Side Rail / Folding Side rail
Bed Side Cabinet, Over Bed Table, and Foot Step
Airmate Cpro Mattress / Urethane

Product Dimension

Dimension
L x W x D : 2230 x 1030 x 760~1005

Product Weight

Weight : ± 90KG
Load Capacity : ± 400KG (Static)
& ± 170KG (Dynamic)
Packing Dimension : 2050 x 1070 mm x 190

Lot No :
Type :

PT. Chitose Internasional Tbk
Jln Industri III No.5 Utama Cimahi
40533 Jawa Barat Indonesia

Showroom - Marketing Office
Distributor Center Jln. H.M.S.
Mina Redja Baros Cimahi 40521
Indonesia Ph +62 22 206 4777

Manufactured by Chitose

245 Kg
175 Kg

MANUFACTURE LABEL AND LOT NUMBER



IZIN EDAR KEMENKES
bahan dasar : sticker warna putih
posisi penempatan pada Frame

Chitose®

MANUFACTURER LABEL



SNI LABEL

