



CB 3003 D-ST

OPTIONS



CB-3003 D ST

FEATURES & PROCEDURES




Head/foot board panel
dismantle

Head/foot board panel
Installation

OPERATION SYSTEM

A unique eye-catching foldable crank is designed to give high-low bed adjustment.




Crank operation system

Side rail compartment
available in every
nursing unit

CASTER OPERATION SYSTEM

Independent lock caster lone at front and one at back of the bed provides stability and secure parking.



PRODUCT SPECIFICATIONS

Operating System Adjustable Hi-Lo, back raise, and foot raise manual by crank

Accessories Side Rail (CSR-001, CSR-002, CSR-003) & Folding Side Rail (FSR-001)

Lock System Double Lock System

Weight (kg) Approx. 90 / **Frame Finishing** Powder Coating

Frame Colour Ivory / **Frame Material** Steel

Dimension (mm) Bed (W 932, L 2141, H 654-994, H/F Board (W 932)

PRODUCT DIMENSION



SNI 12-0179-1987

Type

Lot No.

PT.Chitose Indonesia Mfg.,

Jl. Inggreh No. 8 Bandung - Indonesia
Telp. (022) 6031900 - Fax (022) 6031855

70 mm

Chitose Bed

KEMENKES RI AKD : _____

20 mm

MADE IN INDONESIA

Keterangan :

1. Bahan dasar sticker warna putih
2. Posisi penempelan pada frame